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
**ANALYSIS application laboratory**  
The laboratory is equipped with an ion microscope, an Auger microscope, an Auger/ESCA spectrometer, a microprofiler and an IBS system (Fig. 20).

Fig. 20

### SIMS ANALYSES

Riber offers the following systems:

- the SMS 156 for static SIMS;
- the MIQ 256, under ultra-high vacuum with a quadrupole mass spectrometer;
- a range of combined SIMS/AES/XPS systems of the LAS series.

### ION MICROPROBE MIQ 256

Secondary ion mass spectrometry is a highly sensitive technique for both surface analysis and depth concentration profile determination. It also enables hydrogen to be analysed and can be performed on insulating materials.

To obtain optimum results, the sensitivity must be achieved through:

- clean, ultrahigh vacuum. The MIQ 256 operates at  $7 \times 10^{-10}$  Torr during analysis;
- use of a sensitive mass spectrometer;
- mass filtering the primary ion beam.

Figure 20 shows the sensitivity obtained in a concentration profile of magnesium in silicon. The detection limit is  $4 \times 10^{-11}$  at. cm<sup>-1</sup>.


Fig. 19

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## THE SIA SERIES: AUGER AND ESCA ANALYSIS

The SIA (Surface & Interface Analyzer) systems are optimized for Auger, ESCA and ELS. They are of modular design and include:

- an ultra high vacuum chamber;
- an introduction chamber (several models available);
- a sample manipulator;
- an electron gun for Auger and ELS analysis (several models available);
- an X-ray source for ESCA;
- a computer control system for data acquisition and processing.

## NANOSCAN 50

The Ribet Nanoscan 50 is the only available Auger Microprobe simultaneously featuring a top level electron source, a high performance energy analyzer, and sophisticated computer controlled automation.

It offers:

- very high spatial resolution in SEM/SAM, thanks to a low energy, high brightness, electrostatically focused field emission gun of unique design which delivers a primary current of 1 nA into a 500 Å spot at 15 keV primary energy (10 nA into 1000 Å at 15 keV).
- high sensitivity and chemical information with the high-resolution, high-transmission, MAC 2 electron energy analyzer.
- ease of use and reproducibility through computer control of electron gun, digital scanning unit, electron energy analyzer and motorized sample holder (2" diameter), enabling automatic overnight operation. Deported processors are used to control the system and acquire spectra and SEM/SAM images. This enables true foreground/background operation for the main computer.

Great care has been taken with the system's ergonomics. Two work stations have been designed:

- the SEM/SAM control board commands electron gun parameters, scanning parameters and the digital image acquisition unit. This board, though interfaced to the computer, is manipulated in an analogical manner making system use simple and rapid.
  - the spectra acquisition and data treatment station houses the computer and peripherals. A high level of interaction is achieved through the mouse and powerful software featuring:
    - control of system electronics
    - control of the deported acquisition unit for transfer of spectra to the main computer, and processing
    - transfer of SEM or SAM images acquired by the digital image acquisition unit to the computer, and processing.
- An x-ray source for XPS and ion sources for sputtering are available as options.




Fig. 14

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